


Driva™ Drills Usage Chart

TM, TSV® and Eztetic®


Use a dry erase marker to track the total number of osteotomy preparations that are performed with each drill. It is recommended that reusable drills be replaced after 15 uses. *Single use drills should be disposed of after each procedure.



*Single use drills must be disposed of after each procedure. This chart is intended to be used as a reference guide only. The number of uses per drill will vary and depends on a variety of factors including bone density encountered, proper handling and cleaning. Refer to Instructions for Use (IFU) for Biomet 3i Kits and Instruments (P-ZBDINSTRP) available at ifu.biomet3i.com for complete Cleaning and Sterilization Information.

NOTE: Surgical instruments are susceptible to damage for a variety of reasons including but not limited to: prolonged use, misuse, rough or improper handling, inadequate cleaning, etc. Care must be taken to avoid compromising the intended performance of the instrument.

KNOW HOW TO IDENTIFY EXCESSIVE WEAR



Dull Cutting Edge

Rounded Tip

Rounded Corner

Inspection for Wear and Damage

Before each use, visually inspect all devices and trays/kits for completeness, damage and/or excessive wear (e.g. corrosion or rust build-up on the instrument surface, structural wear or damage, grommet presence or damage partial or complete fracture. If damage or wear is noted that may compromise the function of the instrument or tray, **DO NOT USE**. Contact your Zimmer Biomet representative for replacement.

Refer to Instructions for Use (IFU) for Biomet 3i Kits and Instruments (P-ZBDINSTRP) available at ifu.biomet3i.com for complete Cleaning and Sterilization Information.

