

Screw Fixation System

Quality & Confidence – All In One



1 Secure & Simple

- Easy and stable screw pick-up
- Reliable driver to screw connection
- Stored and sterilized in one durable kit tray

2 Compact & Complete

- Autoclavable metal storage tray
- Heavy, stainless steel screw driver handle
- 1.5 mmD Micro and 2.0 mmD Mini System

3 Versatile Module – Mix & Match

- Titanium alloy fixation screws in 1.5 mmD, 2.0 mmD and 2.3 mmD
- Color coded drills, drivers and screws
- Choose and add for your individual practice needs at any time



ZIMMER BIOMET
Your progress. Our promise.®

Quality & Confidence – All In One

The **Screw Fixation System** provides a compact solution for the temporary fixation and stabilization of bone transplants, suitable resorbable and non-resorbable bone replacement materials, and membranes for the alveolar ridge. Two color coded systems in 1.5 mmD Micro (**Blue**) and 2.0 mmD Mini (**Red**) offer concise and cost effective functional options. The color coding scheme for the two systems, the components and the screws makes rapid identification of the parts possible and simplifies parts matching. This modular storage system permits individual configuration and the open design ensures optimum access during cleaning and sterilization. Fixation screws and mesh are manufactured from pure titanium or titanium alloy. They are biocompatible, corrosion-proof and nontoxic in the biological environment. They allow imaging virtually free of artifacts.



Secure & Simple

The screws are easily picked up and have a reliable connection to the driver.



Getting Started Made Easy With Our Start-Up Screw Fixation Kit



Start-Up Screw Fixation Kit* (69.01.10Z) features:

- Durable, autoclavable metal storage tray
- Heavy stainless steel screw driver handle
- Micro System, 1.5 mmD
 - Stainless steel screw driver inserts (long & short)
 - Stainless steel pilot drill, 14 mmL
 - Stainless steel pilot block drill, 15 mmL

Screw Fixation kit – Start-Up Kit Ordering Information

Contents	Tray	Screw Driver Handle	Screw Driver Insert, Short	Screw Driver Insert, Long	Pilot Drill, Micro, 14 mmL	Pilot Block Drill, Micro
Assembled Start-Up Kit*	69.01.10Z					
Components	69.01.11Z	75.23.52Z	75.23.23Z	75.23.19Z	69.01.09Z	69.01.16Z

*Screws are not included in the start-up kit

Contact us at 1-800-342-5454 or visit
zimmerbiometdental.com

Zimmer Biomet Dental
 Global Headquarters
 4555 Riverside Drive
 Palm Beach Gardens, FL 33410
 Tel: +1-561-776-6700
 Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. Screw Fixation System is manufactured by Medicon e.G. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0205 REV B 09/19 ©2019 Zimmer Biomet. All rights reserved.

