

RealGUIDE™ Z3D Guided Surgery Kit

Designed for Zimmer Biomet's TSV®, Trabecular Metal™ and Eztetic® Dental Implants

FULLY
GUIDED!



 RealGUIDE™ Z3D
GUIDED SURGERY KIT



RealGUIDE Z3D GUIDED SURGERY KIT

The RealGUIDE Z3D Guided Surgery Kit contains instruments needed for the planned surgery. It is designed to work with Zimmer Biomet's **TSV, Trabecular Metal and Eztec** Implant Systems.

The kit information is available in the RealGUIDE Software Suite, an open system for diagnosis and advanced surgical guide design that will support you with customized drilling reports to be well-prepared for dental implant therapy.

- Fully guided placement of Zimmer Biomet Dental Eztec 3.1 mm Implant and TSV/TM Implants in sizes 3.7 mm, 4.1 mm and 4.7 mm
- Color-coded for ease of use
- Compact kit, which is useful for space-constrained dental offices
- No handles or keys are necessary while preparing the osteotomy
- DLC (Diamond-like Carbon) coating on the drills provides wear resistance characteristics
- Includes implant mounts hex feature to align with the hex feature on the corresponding guide sleeves
- Staging pedestal for holding the implant vial to remove the Fixture Mount Transfer (FMT) of the TSV implants is available. This enables the TSV implants to be placed fully guided through the surgical guides
- Both soft and dense bone protocols are available



RealGUIDE Software Suite
with PLAN, APP, GUIDE Modules

Product	Part No.
RealGUIDE Z3D Surgical Kit	3DM0070Z



MORE INFORMATION
ABOUT THE RealGUIDE Z3D
GUIDED SURGERY KIT AT
ZBDGUIDEDSURGERY.COM

Contact us

at 1-800-342-5454 or visit zbdguidedurgery.com

 **RealGUIDE™ Z3D**
GUIDED SURGERY KIT

All content herein is protected by copyright, trademarks, and other intellectual property rights, as applicable, owned by or licensed to 3Diemme SRL or its affiliates unless otherwise indicated, and must not be redistributed, duplicated, republished or reprinted, in whole or in part, without express written consent of the owner. Unless otherwise indicated, all products are manufactured by 3Diemme SRL or its affiliates, and marketed and distributed by 3Diemme SRL and its authorized marketing partners. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries and/or regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. Distribution to any other recipient is prohibited. ZB1317 REV A 04/21 ©2021 Zimmer Biomet. All rights reserved.

CE
0051 Certificate n.003/MDR

