

# Prosthodontic Screwdriver



Manufacturer:

W&H Dentalwerk Bürmoos GmbH

Distributor:

**BIOMET 3i™**  
PROVIDING SOLUTIONS - ONE PATIENT AT A TIME™

# Prosthodontic Screwdriver

The cordless Prosthodontic Screwdriver has torque control, which enables precise tightening of retaining screws. This cordless screwdriver reduces the need for traditional mechanical rotary ratchet wrench systems.



## Features Include:

- Torque control of 8-40 Ncm in 1 Ncm steps for precise tightening of retaining screws
- An automatic stop when the preset torque is reached
- Cordless operation with constant drive speed
- A small contra-angle head for better visibility of the treatment site
- A contra-angle that can be rotated 360°
- A clockwise and counter clockwise operation
- A sterilizable contra-angle

## Prosthodontic Screwdriver IA-40 Consists Of:

- A handpiece drive IA-40H with keypad and rechargeable Li-Ion battery (not sterilizable)
- A contra-angle IA-80 with a hexagon chucking system
- A charging station with country specific cables available

## Technical Data

Handpiece Drive IA-40H	Contra-Angle Handpiece IA-80
Rechargeable Battery Type: Li-Ion	80:1 with W&H Special Coupling
Battery Voltage: 3.7 V	Weight: 55 g/1.9 oz
Nominal Capacity: 680 mAh	<b>Charging Station</b>
Speed on the Instrument: 25 rpm	Supply Voltage: 100–240 V
Torque: 8–40 Ncm	Frequency: 50–60 Hz
Weight: 90 g/3.2 oz	Rated Current: 0.08–0.12 A
Charging Time: Approx. 100 min	Power: 5 W
Battery Capacity: Approx. 40 screws at medium torque	Weight: 345 g/12.2 oz

Distributor:

**BIOMET 3i**

Check Out The **BIOMET 3i** Online Store at [www.shopbiomet3i.com](http://www.shopbiomet3i.com)

(Available for Spain & US Customers Only)



Global Headquarters  
4555 Riverside Drive  
Palm Beach Gardens, FL 33410  
Toll Free: 1-800-342-5454  
Phone: +1-561-776-6700  
Fax: +1-561-776-1272  
[www.biomet3i.com](http://www.biomet3i.com)

Dental Ibérica S.L.  
WTC Alameda Park, Ed. 4, Planta 2  
C/Tirso de Molina, 40  
08940, Cornellà de Llobregat  
(Barcelona) Spain  
Phone: +34-93-470-55-00  
Fax: +34-93-371-78-49  
[www.biomet3i.com](http://www.biomet3i.com)

**Manufacturer:**  
W&H Dentalwerk Bürmoos GmbH  
Ignaz-Glaser-Straße 53  
5111 Bürmoos, Austria  
Phone: +43-6274-6236-0  
Fax: +43-6274-6236-55  
[www.wh.com](http://www.wh.com)



BIOMET 3i is a trademark of BIOMET 3i LLC. W&H is a registered trademark of W&H Dentalwerk Bürmoos GmbH.  
©2015 BIOMET 3i LLC.

This material is intended for clinicians only and is NOT intended for patient distribution. This material is not to be redistributed, duplicated or disclosed without the express written consent of BIOMET 3i.



ART1241  
REV A 06/15