

Puros[®] DBM Particulate



1 Osteoinductive potential preserved¹

- Preservation of protein activity maintains osteoinductivity (OI) potential^{2,3}

2 Best in class processing¹

- Proprietary Cancellé SP[®] processing designed to preserve protein activity allowing the DBM to maintain OI potential
- Only DBM with verified OI potential is released to processing
- Low temperature, low dose gamma irradiation preserves the utility of the graft

3 Trusted Source¹

- Safety you expect from the Puros Brand
- Single donor source
- Terminally sterilized to SAL 10⁻⁶

4 Convenient Handling

- Quick hydration
- Room temperature storage⁴
- 2 year shelf-life

Ordering Information:

Description	Volume	Item No.
Puros DBM Particulate 0.25 cc	0.25 cc	DBMPART025
Puros DBM Particulate 0.5 cc	0.5 cc	DBMPART050
Puros DBM Particulate 1.0 cc	1.0 cc	DBMPART100

1. Data on file with RTI Surgical.
2. Each lot of Puros DBM particulate is tested for OI potential using an in vivo athymic rat model.
3. Performance data from animal studies may not be representative of performance in humans.
4. Refer to implant labeling for specific storage conditions.



Contact us at 1-800-342-5454 or visit

zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Tel.: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. Cancellé SP is a U.S. registered trademark of RTI Surgical, Inc. Puros Allografts are manufactured by RTI Surgical, Inc. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB1195 REV A 11/20 ©2020 Zimmer Biomet. All rights reserved.

