

RegenerOss® Allograft Putty Plus

Graft With Ease!



1 Ideal Handling

The lecithin carrier provides a moldable vehicle for containing and transporting allograft bone granules.

2 Osseinductive & Osseoconductive

RegenerOss Allograft Putty Plus is comprised of a mixture of DBM and mineralized allograft to facilitate the regeneration of new bone.

3 Simple Syringe Delivery

The ergonomic design features a small diameter syringe with a curved tip for easy access to defect sites.



ZIMMER BIOMET
Your progress. Our promise.™

RegenerOss® Allograft Putty Plus

ENCOURAGES BONE GROWTH POTENTIAL BY INCORPORATING A HIGH BONE-TO-CARRIER RATIO WITHOUT SACRIFICING ITS HANDLING CHARACTERISTICS.¹

CARRIER

IDEAL HANDLING

- Resists movement during irrigation and keeps bone particles from migrating
- Easy to mold – doesn't break apart
- Potentially enhances osseoinductive properties of the Demineralized Bone Matrix (DBM)²

BONE

OSSEOINDUCTIVE AND OSSEOCONDUCTIVE

- Osseoinductive potential of DBM encourages bone growth¹
- Mineralized cancellous chips provide a natural scaffold to maintain space and allow for attachment of cells
- 48% bone graft material by weight
- Demineralized PLUS Mineralized Bone³
- 28% Cortical Demineralized Bone Matrix
- 20% Cancellous Mineralized Bone Chips

SIMPLE SYRINGE DELIVERY

- Enables clinicians to treat hard-to-reach defects
- The plunger extends all the way to the end of the syringe for complete expulsion of the putty, minimizing waste


Additional Benefit:

Mineralized component is radiopaque for post-operative verification of grafting

Indications:

- Peri-implant defects
- Localized ridge augmentation
- Extraction sockets
- Periodontal defects
- Sinus augmentation
- Cystic defects
- Defects associated with root resection or apicoectomy

RegenerOss
Putty Plus



Ordering Information:

Item Number	Description
ROAPM05	05. cc
ROAPM10	1.0 cc
ROAPM20	2.0 cc (Four 0.5 cc Syringes)

Value
Pack!

1. Urist, MR. Bone: Formation By Autoinduction. Science, 1965, 150: pp 893.
2. Han, B, et al. Combined effects of phosphatidylcholine and demineralized bone matrix on bone induction. Connective Tissue Research, 2003, 44: pp 160–166.
3. Tissue manufactured and provided by LifeLink Tissue Bank, an AATB certified tissue bank.

Contact us at 1-800-342-5454 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Phone: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc., and distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). RegenerOss Allograft is manufactured by Interpore Cross, Inc. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0061 REV A 02/17 ©2017 Zimmer Biomet. All rights reserved.

