

# Screw Fixation System

## Quality & Confidence – All In One



The **Screw Fixation System** provides a compact solution for the temporary fixation and stabilization of bone transplants, suitable resorbable and non-resorbable bone replacement materials, and membranes for the alveolar ridge. Two color coded systems, 1.5 mmD MICRO (blue) and 2.0 mmD MINI (red), make easy and rapid identification of the parts possible and simplifies parts matching.

This modular storage system permits individual configuration and the open design aims to ensure convenient access during cleaning and sterilization.

Fixation screws and mesh are manufactured from pure titanium or titanium alloy. They are biocompatible, corrosion-proof and nontoxic in the biological environment. They allow imaging virtually free of artifacts.



### Secure & Simple

The **Screw Fixation System** comes with power grip for secure and stable transfer to the surgical site. The screws are easily picked up and have a reliable connection to the driver.

### Well-Arranged

The screw cartridges of the **Screw Fixation System** are more than simply storage containers. The cartridges make organization and sterilization of screws easy. An imprinted four data identification system provides traceability.

# Screw Fixation Kit – Start-Up Kit Ordering Information




Contents	Tray	Screw Driver Handle	Screw Driver Insert, Short	Screw Driver Insert, Long	Pilot Drill, Micro, 14 mmL	Pilot Block Drill, Micro
Assembled Start-Up Kit*	69.01.10Z					
Components	69.01.11Z	75.23.52Z	75.23.23Z	75.23.19Z	69.01.09Z	69.01.16Z

\*Screws are not included in the start-up kit



1.5 mmD MICRO System	System/Color Indication	Color Marking	Measurements	Kit & Stand Alone	Stand Alone Only
Tray	1.5 Micro & 2.0 Mini; Blue & Red	n/a	n/a	69.01.11Z	
Screw Driver Handle	1.5 Micro & 2.0 Mini; Blue & Red	n/a	n/a	75.23.52Z	
Screw Driver Insert, Short	1.5 Micro; Blue	blue	17 mmL	75.23.23Z	
Screw Driver Insert, Long	1.5 Micro; Blue	blue	55 mmL	75.23.19Z	
Pilot Drill, Micro, Short, 9 mmL	1.5 Micro; Blue	blue	1.2 mmD/9 mmL		69.01.14Z
Pilot Drill, Micro, Long, 14 mmL	1.5 Micro; Blue	blue	1.2 mmD/14 mmL	69.01.09Z	
Pilot Block Drill, Micro	1.5 Micro	green	1.6 mmD x 15 mmL	69.01.16Z	




Fixation Screws 1.5 mmD MICRO System	Cartridge (10 screws)	System	Screw Color	Measurements Length mm	Micro 1.5 mmD Self drilling	Micro 1.5 mmD Self tapping
Screw, self drilling	1	1.5 Micro	green	3.5	68.85.83Z	---
Screw, self drilling	1	1.5 Micro	green	4	68.85.84Z	---
Screw, self drilling	1	1.5 Micro	green	5	68.85.85Z	---
Screw, self drilling	1	1.5 Micro	green	7	68.85.87Z	---
Screw, self tapping	1	1.5 Micro	silver	9	---	68.85.49Z
Screw, self tapping	1	1.5 Micro	silver	11	---	68.85.51Z
Screw, self tapping	1	1.5 Micro	silver	13	---	68.85.53Z
Screw, self tapping	1	1.5 Micro	silver	15	---	68.85.55Z
Screw, self tapping	1	1.5 Micro	silver	17	---	68.85.57Z



2.0 mm MINI System	System/Color Indication	Color Marking	Measurements	Kit & Stand Alone	Stand Alone Only
Tray	1.5 Micro & 2.0 Mini; Blue & Red	n/a	n/a	69.01.11Z	
Screw Driver Handle	1.5 Micro & 2.0 Mini; Blue & Red	n/a	n/a	75.23.52Z	
Screw Driver Insert, Short	2.0 Mini; Red	red	17 mmL		75.23.21Z
Screw Driver Insert, Long	2.0 Mini; Red	red	55 mmL		75.23.22Z
Pilot Drill, Mini, Short	2.0 Mini; Red	red	1.6 mmD/9 mmL		69.01.15Z
Pilot Block Drill, Mini	2.0 Mini	yellow	2.1 mmD x 15 mmL		69.01.17Z



Fixation Screws 2.0 mm MINI System	Cartridge (10 screws)	System/Color Indication	Screw Color	Measurements Length mm	Micro 2.0 mmD Self drilling
Screw, self tapping	1	2.0 Mini; Red	silver	4	68.94.00Z
Screw, self tapping	1	2.0 Mini; Red	silver	5.5	68.94.01Z
Screw, self tapping	1	2.0 Mini; Red	silver	7	68.94.02Z
Screw, self tapping	1	2.0 Mini; Red	silver	9	68.94.03Z
Screw, self tapping	1	2.0 Mini; Red	silver	11	68.94.04Z
Screw, self tapping	1	2.0 Mini; Red	silver	13	68.94.05Z
Screw, self tapping	1	2.0 Mini; Red	silver	15	68.94.06Z
Screw, self tapping	1	2.0 Mini; Red	silver	17	68.94.07Z



Additional Components	Color Indication	Color Marking	Measurements	Stand Alone Only	
Titanium Mesh	Blue & Red	n/a	42x32 mm; 0.1 mm thick	68.86.95Z	
Titanium Mesh	Blue & Red	n/a	42x32 mm; 0.2 mm thick	68.86.96Z	
Screw, Rescue	1 Cartridge (10 screws)	Red	pink	7 mm	68.94.10Z
Screw, Rescue	1 Cartridge (10 screws)	Red	pink	11 mm	68.94.11Z

Contact us at 1-800-342-5454 or visit  
[zimmerbiometdental.com](http://zimmerbiometdental.com)

Zimmer Biomet Dental  
 Global Headquarters  
 4555 Riverside Drive  
 Palm Beach Gardens, FL 33410  
 Tel: +1-561-776-6700  
 Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc. and marketed and distributed by Zimmer Biomet Dental and its authorized marketing partners. Screw Fixation System is manufactured by Medicon e.G. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations. This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental. ZB0204 REV B 09/19 ©2019 Zimmer Biomet. All rights reserved.

