

Puros[®] Cortico-Cancellous Particulate Allograft



The Best of Both Worlds.

1 The Natural Mix For Healthy Bone Growth

Puros Cortico-Cancellous Particulate Allograft has been designed to provide an anatomic-based mix of 70% cortical and 30% cancellous bone particulate to simplify your grafting process

3 A Trustworthy Source

Safety and performance you expect from the Puros brand

Single donor vials

Two particulate sizes to accommodate different defect sizes

2 Easy To Use

Retains osteoconductive properties due to the preservation of the natural bone matrix collagen and mineral composition, trabecular pattern and original porosity¹


Easy handling—quick hydration, five-year shelf life and room-temperature storage

No need to mix on site

Puros Cortico-Cancellous Particulate Allograft

Sterilized and preserved using the proprietary Tutoplast® process, Puros Cortico-Cancellous Particulate Allograft is a high-quality allograft designed for large and small volume bone regeneration procedures.

Accuracy Assured Convenient Packaging




The Unique Tutoplast Process


The proprietary Tutoplast process assures a high standard of tissue safety and quality with minimal risk of disease transmission.

The process preserves the valuable collagen matrix and tissue integrity while inactivating pathogens and gently removing unwanted materials, such as cells, antigens and viruses. The result is safe, biocompatible tissue.


For over 40 years, a variety of Tutoplast processed tissues have been safely used in more than five million procedures.¹




Delipidization




Osmotic treatment



Oxidative treatment



Solvent dehydration



Low-dose gamma irradiation

Ordering Information

Catalog Number	Description
68800	Puros Cortico-Cancellous Particulate, 0.5 cc, 250-1000 µm
68801	Puros Cortico-Cancellous Particulate, 1 cc, 250-1000 µm
68802	Puros Cortico-Cancellous Particulate, 2 cc, 250-1000 µm
68803	Puros Cortico-Cancellous Particulate, 0.5 cc, 1000-2000 µm
68804	Puros Cortico-Cancellous Particulate, 1 cc, 1000-2000 µm
68805	Puros Cortico-Cancellous Particulate, 2 cc, 1000-2000 µm

Zimmer Biomet Dental offers a comprehensive line of allografts for bone augmentation needs.

¹ Data on file with RTI Surgical, Inc.

Contact us at 1.800.363.1980 or visit zimmerbiometdental.com

Zimmer Biomet Dental
Global Headquarters
4555 Riverside Drive
Palm Beach Gardens, FL 33410
Tel: +1-561-776-6700
Fax: +1-561-776-1272

Unless otherwise indicated, as referenced herein, all trademarks are the property of Zimmer Biomet; and all products are manufactured by one or more of the dental subsidiaries of Zimmer Biomet Holdings, Inc., and distributed and marketed by Zimmer Biomet Dental (and, in the case of distribution and marketing, its authorized marketing partners). Tutoplast is a U.S. registered trademark owned by Tutogen Medical GmbH. Puros products are manufactured by RTI Surgical, Inc. dba RTI Biologics or Tutogen Medical GmbH. For additional product information, please refer to the individual product labeling or instructions for use. Product clearance and availability may be limited to certain countries/regions. This material is intended for clinicians only and does not comprise medical advice or recommendations.

This material may not be copied or reprinted without the express written consent of Zimmer Biomet Dental.

ZB0018CAN REV B 03/18 ©2018 Zimmer Biomet. All rights reserved.

